

Equine 90/30 Portable Radiographic Technique Chart

STUDY	VIEW	SID	kVp	TIME	mA	mAs
NAVICULAR	65 Degree Navicular	30"	80	0.06	20	1.20
	65 Degree P-3	30"	80	0.06	20	1.20
	DP	30"	80	0.06	20	1.20
	Lateral	30"	80	0.06	20	1.20
	Flexor Skyline	30"	80	0.06	20	1.20
	Navicular Obliques	30"	80	0.06	20	1.20
	65 P-3 Obliques	30"	80	0.03	20	0.60
P3	Lateral	30"	80	0.06	20	1.20
	65 Degree P-3	30"	80	0.03	20	0.60
	65 P-3 Obliques	30"	80	0.03	20	0.60
PASTERN	DP	30"	80	0.06	20	1.20
	Lateral	30"	80	0.06	20	1.20
	Medial Oblique	30"	80	0.06	20	1.20
	Lateral Oblique	30"	80	0.06	20	1.20
FETLOCK	DP	30"	80	0.06	20	1.20
	Medial Oblique	30"	80	0.06	20	1.20
	Lateral	30"	80	0.06	20	1.20
	Lateral Oblique	30"	80	0.06	20	1.20
	Flexed Lateral	30"	80	0.06	20	1.20
	Special Palmer	30"	90	0.06	15	0.90
METACARPUS/METATARSUS	DP	30"	80	0.06	20	1.20
	Lateral	30"	80	0.06	20	1.20
	Medial Oblique	30"	80	0.06	20	1.20
	Lateral Oblique	30"	80	0.06	20	1.20
CARPUS	DP	30"	80	0.06	20	1.20
	Medial Oblique	30"	80	0.06	20	1.20
	Lateral	30"	80	0.06	20	1.20
	Lateral Oblique	30"	80	0.06	20	1.20
	Flexed Lateral	30"	80	0.06	20	1.20
	Distal Rad. Sky	30"	80	0.06	20	1.20
	Distal Row Sky	30"	80	0.06	20	1.20
	Proximal Row Sky	30"	80	0.06	20	1.20

STUDY	VIEW	SID	kVp	TIME	mA	mAs
TARSUS	DP	30"	80	0.06	20	1.20
	Medial Oblique	30"	80	0.06	20	1.20
	Lateral	30"	80	0.06	20	1.20
	Lateral Oblique	30"	80	0.06	20	1.20
	Calcaneal Skyline	30"	80	0.06	20	1.20
SKULL TMJ & BULLAE	DV	30"	90	0.25	9	2.25
	Lateral	30"	90	0.15	9	1.35
	Right Oblique	30"	90	0.15	9	1.35
	Left Oblique	30"	90	0.15	9	1.35
TEETH	All Views	30"	84	0.36	9	3.24
SKULL SINUS & DENTAL	DV	30"	90	0.36	9	3.24
	Lateral	30"	80	0.06	20	1.20
	Right Oblique	30"	80	0.06	20	1.20
	Left Oblique	30"	80	0.06	20	1.20
CERVICAL SPINE	C1, C2/C2,C3	30"	80	0.26	12	3.12
	C4, C5	30"	84	0.30	9	2.70
	C5, C6	30"	86	0.32	9	2.88
	C6, C7	30"	86	0.36	9	3.24
PHARYNX	Lateral	30"	80	0.06	20	1.20
SHOULDER	Lateral	30"	86	0.32	9	2.88
ELBOW	Cranial Caudal	30"	84	0.26	9	2.34
	Medio-Lateral	30"	84	0.30	9	2.70
RADIUS	Cranial Caudal	30"	82	0.20	9	1.80
	Lateral	30"	82	0.20	9	1.80
STIFLE	Caudal Cranial	30"	84	0.30	9	2.70
	Lateral	30"	80	0.26	12	3.12
	Lateral Oblique	30"	80	0.16	12	1.92
	Patellar Skyline	30"	80	0.26	12	3.12
TIBIA	Caudal Cranial	30"	82	0.20	9	1.80
	Lateral	30"	82	0.20	9	1.80
WITHERS	Lateral	30"	80	0.26	12	3.12